



## PRESS RELEASE

18th May, 2009

# NRS II™ Transwab®, new generation of environmental sampling for food and pharma

NRS II™ Transwab® is a new range of pre-filled sampling devices with a choice of fill volume and diluent formulation.

Each device comes as a pre-filled self-standing tube, complete with rayon swab, ready to use in a whole variety of sampling situations. The swab is already mounted in the screw-cap. There are no peel pouches or plugs to dispose of, reducing packaging, and the risk of components getting lost on food or sterile product lines.

The tubes are available with a choice of 1ml, 5ml and 10ml precision fills. There is also a complete range of diluents including NRS™ medium, Buffered Peptone Water, Butterfield Buffer, Dey-Engley Neutralising Broth, Lethen Broth and MRD (also called Peptone Saline).

All NRS II™ Transwab® devices feature Medical Wire's new leak-proof self-standing screw cap tubes.

The advertisement graphic features the mwe medical wire logo at the top left. Below it, the product name 'NRS II™ Transwab®' is displayed in a blue banner. The main text reads: 'A new generation of environmental sampling devices for the food and pharmaceutical industries.' A bulleted list of features includes: 'Pre-filled self-standing tubes', 'New rayon bud retains integrity in solution', 'Optimum sampling, maximum recovery', 'Choice of diluent', and 'Conform to national and international standards'. An image shows a blue screw cap tube with a rayon swab attached. A small inset image shows orange rod-shaped bacteria on a red surface, labeled 'Coliform bacteria'. A vertical blue line with a white bulbous end hangs from the top right. The bottom of the graphic has a blue banner with the text 'A new generation of environmental sampling'.

**mwe**  
medical wire

### NRS II™ Transwab®

A new generation of environmental sampling devices for the food and pharmaceutical industries.

- Pre-filled self-standing tubes
- New rayon bud retains integrity in solution
- Optimum sampling, maximum recovery
- Choice of diluent
- Conform to national and international standards

All NRS II™ Transwab® devices feature Medical Wire's new leak proof labelled self-standing screw cap tube made from shatter proof polypropylene, with a highly visible blue shaft swab attached to the cap.

Coliform bacteria

A new generation of environmental sampling

[www.mwe.co.uk](http://www.mwe.co.uk)  
[www.mwe-usa.com](http://www.mwe-usa.com)

Medical Wire & Equipment, Corsham, Wiltshire, SN13 9RT, England