

In today's fast-paced world, healthcare professionals are constantly seeking efficient solutions to improve patient care and streamline processes. Our self-collection range offers a compassionate and practical approach to specimen collection, making it easier for patients while supporting healthcare providers in their mission to deliver exceptional care.

When it comes to healthcare, there's no place like home. We know that, just like you, your patients want to ensure their families are well taken care of. Our self-collection solutions allow patients to feel more comfortable and empowered by providing privacy and reducing the burden of in-person visits. These options allow for greater flexibility and peace of mind for individuals and families, while also saving valuable time for healthcare professionals. It's about making healthcare simpler, more personal, and, most importantly, more focused on results that matter.

Our self-collection products offer a seamless, reliable experience for both healthcare providers and patients.

Choose the solution that's right for you:
GOOD, BETTER, OR BEST

SELF-COLLECTION

U-SWAB™



Convenience: Collection at the patient's convenience, without compromising quality or accuracy.



Privacy: Patients can collect samples discreetly, empowering them to take charge of their health.



Time & Money-Saving: Efficient processes that reduce the need for time-consuming in-person visits.



Experience a New Level of Comfort: No swabs needed with our innovative SampleSUITE™

3 FORMATS


GOOD, BETTER & BEST



U-SWAB™



U-SWAB™
Liquid Transport



SampleSUITE™
+ U-SWAB™ Liquid Transport

U-SWAB™

SELF-COLLECTION

WHAT YOU NEED TO KNOW



- Consumer instructions for use available on U-SWAB.CO.UK
- Tamper Proof Seal
- Shelf-Life: 24 months

Specifications


U-SWAB™ Liquid Transport Media

Swab Length: 150mm
 Swab Breakpoint: 79mm
 Small Scoop Length: 79mm
 Large Scoop Length: 99mm

U-SWAB™

Swab Length: 120mm
 Swab Breakpoint: 80mm
 Scannable & Peelable Barcode Option Available

| | Range | Product Code | Description |
|---------------|--------------------------------|--------------|--|
| BEST | SampleSUITE™ | M200121 | StoolPOD™ – Blister Pouch |
| | SampleSUITE™ | M200122 | CapillaryPOD™ – 40µl – Blister Pouch |
| | SampleSUITE™ | M200123 | SalivaPOD™ – ECO – Blister Pouch |
| BETTER | U-SWAB™ Liquid Transport Media | M200141 | U-SWAB™ MRSA Trio – 1 Standard Plastic Foam Printed Breakpoint + 1 Blue Plastic Foam + 1 Red Plastic Foam – 12x80 Vial – 1.5ml Liquid Amies Medium |
| | U-SWAB™ Liquid Transport Media | M200142 | U-SWAB™ MRSA Duo -1 Standard Plastic Foam Printed Breakpoint + 1 Red Plastic Foam – 12x80 Vial – 1.5ml Liquid Amies Medium |
| | U-SWAB™ Liquid Transport Media | M200145 | U-SWAB™ Fecal – Scoop – 12x80 Vial – 2ml Cary Blair Medium |
| | U-SWAB™ Liquid Transport Media | M200146 | U-SWAB™ Fecal – Scoop – 16x100 Vial - 2ml Cary Blair Liquid Medium |
| GOOD | U-SWAB™ | M200147 | U-SWAB™ Fecal – Scoop – 12x80 Vial – Dry |
| | U-SWAB™ | M200148 | U-SWAB™ Fecal – Scoop – 16x100 Vial – Dry |
| | U-SWAB™ | M200101 | U-SWAB™ Vaginal – Standard Flockswab – Labelled Tube |
| | U-SWAB™ | M200102 | U-SWAB™ Vaginal – Standard Flockswab – Labelled Tube – Scannable Peel Barcode |



We ensure every sample is collected and transported in optimal condition, delivering results that advance care, discovery, and progress.

HAVE YOU THOUGHT ABOUT:

Giving your patients the flexibility to choose the self-collection device that suits them best – from the simple **U-SWAB™ Vaginal** to the innovative, swab-free **SampleSUITE™** option in our Best tier. Empower patients with personalised care options.

✉ info@mwe.co.uk

🌐 U-SWAB.CO.UK
mwe.co.uk

🌐 Follow us on LinkedIn